

## **FLEXILOG pH Oesophageal pH recorder**



- **Small: 12.5 x 9 x 3cm (L x W x D)**
- **Lightweight: 250g including batteries**
- **Large screen LCD**
- **Up to 4 channels of pH**
- **Optional pressure channel for LES detection**
- **4 patient event buttons**
- **Comfortable carrying case**
- **Minimum study preparation time**
- **All actions prompted by LCD**
- **Compatible with antimony and glass pH catheters**
- **LES location with water perfused or solid state transducer**
- **USB communication to computer running FLEXISOFT III**

Oakfield Instruments Ltd  
Oakfield Industrial Estate  
Eynsham, Witney  
Oxon, OX29 4AN, UK

T: +44 (0) 1865 882532  
F: +44 (0) 1865 883970  
E: [info@flexilog.net](mailto:info@flexilog.net)