

FLEXILOG

Gastrointestinal Motility Products

Oesophageal pH Catheters



FLEXILOG

Built for Reliability and accuracy, FLEXILOG catheters are suitable for use with any FLEXILOG or other ambulatory pH System.

Disposable Catheters (Internal Reference);

- Suitable for adult or pediatric use
- 1.8mm diameter
- Marked every cm with a dot and numbered every 5cm
- Single Use
- Pack of 10

Single channel		14713
Dual Channel	5cm spacing	14714
	10cm spacing	14715
	15cm spacing	14716

Reusable Catheters (Internal Reference);

- 2.1mm diameter
- Marked every cm with a dot and numbered every 5cm
- Pack of 5

Single channel		14246
Dual Channel	5cm spacing	14247
	10cm spacing	14248
	15cm spacing	14249

Reusable Catheters (External Reference);

- 2.1mm diameter
- External skin reference
- Marked every cm with a dot and numbered every 5cm
- Pack of 5
- Includes 10 adhesive rings

Single channel		14531
Dual Channel	3cm spacing	14532
	5cm spacing	14533
	7cm spacing	14534
	10cm spacing	14535
	15cm spacing	14536

Paediatric Reusable Catheters (External Reference);

- **Very thin 1.5mm diameter**
- External Skin Reference
- Marked every cm with a dot and numbered every 5cm
- Pack of 5
- Includes 10 adhesive rings

Single channel		14527
Dual Channel	3cm spacing	14528
	5cm spacing	14529
	7cm spacing	14530



Buffers

- For use with any FLEXILOG pH Catheters
- 250ml Bottle

pH 1.07	13847
pH 4.00	13848
pH 7.01	13849



External Skin Reference Accessories

- For use with any FLEXILOG External Reference pH Catheter

Adhesive Rings (Pack of 10)	14558
Electro-Conductive gel	14559

Oakfield Instruments Ltd
Oakfield Industrial Estate
Eynsham, Witney
Oxon OX29 4AN
United Kingdom

T +44 (0)1865 882532
F +44 (0)1865 883970
E info@flexilog.net



www.flexilog.net

Oakfield Instruments Inc
8917 Glenwood Avenue
Raleigh NC 27617
USA

T +1 919 571 0110
F +1 775 305 1431
E info@flexilog.us