


Disposable Bipolar Hemostasis Probe



US Endoscopy offers a new twist in bipolar hemostasis probes. Our 7FR and 10FR probes contain a hydrophilic coating called Hydromer® that is grafted to the tips of these devices to help prevent tissue from sticking to the probes. Their extra long, 350cm length allows breathing room between the physician and nurse, and flexibility in the arrangement between the electrosurgical generator, patient and endoscopy team. US Endoscopy bipolar probes are ideal for use in EGD, enteroscopy and colonoscopy procedures – simplifying materials management.



- **Universal 350cm length**

- *fits EGD scopes, enteroscopes and colonoscopes*
- *allows breathing room between physician and nurse*
- *eliminates need to reposition the generator, patient or endoscopy team as is often necessary with shorter bipolar probes*

- **Hydromer® grafted tip**

helps prevent tissue from sticking to probe

- **Various sizes and electrosurgical connector options**

- **Realistically priced**





Disposable Bipolar Hemostasis Probe



ORDERING INFORMATION 800-769-8226

| Product Number | Description | Size | Length (cm) |
|----------------|------------------------------------|-------|-------------|
| 00711841 | Valley Lab, Everest Medical BiCoag | 7 FR | 350 |
| 00711843 | BICAP, BICAP II, Bard, Microvasive | 7 FR | 350 |
| 00711847 | BICAP, BICAP II, Bard, Microvasive | 10 FR | 350 |

Packaged as 1 each



5976 Heisley Road, Mentor, Ohio 44060
Phone (440) 639-4494 • Fax (440) 639-4495
Customer Service (800) 769-8226
Email: info@usendoscopy.com www.usendoscopy.com