



134 °C
273 °F

Mechanical Lithotripter Set Type III

- Combination of methodical and emergency lithotripter in one device.
- Easy handling by unique spiral and traction cable fixation.
- Compatible to all common ENDO-FLEX Lithotripsy baskets.



Order No.

Mechanical Lithotripter Set Type III

120 610	1 x storage case 1 x Lithotripsy handle 2 x stainless steel spiral Ø 2.6 mm / length 195 cm 1 x stainless steel spiral Ø 4.0 mm / length 100 cm 1 x basket, Ø 2.6 mm, 4 wires, basket type, size 50 mm 1 x basket, Ø 2.6 mm, 4 wires, basket type, size 60 mm 1 x basket, Ø 2.6 mm, 4 wires, basket type, size 70 mm
---------	--

Alternatively to above complete set also available:

Order No.

120 611	Lithotripsy handle 1 x Spiral Ø 2.6 mm / length 195 cm 1 x Spiral Ø 4.0 mm / length 100 cm
---------	--

reusable

Mechanical Lithotripter Set

134 °C
273 °F

Mechanical Lithotripter Set Type I

- for methodical lithotripsy "Through-The-Scope"



Order No.	Mechanical Lithotripter Set Type I
120 265N	1 x handle with coaxial fixation 1 x stainless steel spiral Ø 2.6 mm 2 baskets, Ø 2, 6, 4 wires, basket type, size 50 mm
120 266N	1 x handle with coaxial fixation 1 x stainless steel spiral Ø 2.6 mm 2 baskets, Ø 2, 6, 4 wires, basket type, size 60 mm
120 267N	1 x handle with coaxial fixation 1 x stainless steel spiral Ø 2.6 mm 2 baskets, Ø 2, 6, 4 wires, basket type, size 70 mm