

PlumePort™ SEO

Laparoscopic Plume Filtration Device

Features:

Controlled Management

of plume throughout the
surgical procedure

Universal Connection

Standard luer lock fitting and
connects to most surgical
plume evacuators

Approximately 50% More Flow

than passive laparoscopic filters while
maintaining active suction and
peritoneal distention

Quiet Operation

Adjustable Flow Rate

Via surgical plume evacuator

Excellent Value

An economical effective solution

Easy to Use

Latex free sterile tubing



PlumePort™ SEO easily connects to standard trocars to safely remove the surgical plume during laparoscopy. Utilizing available surgical plume evacuators and coupling it with this calibrated device, the PlumePort SEO provides fast and safe plume management while maintaining peritoneal distention without OR personnel interaction.

www.buffalofilter.com

Technical Information

Flow Rate

ViroVac on Medium	(8.3 lpm)
PlumeSafe™ TURBO (60%)	(5.1 lpm)
PlumeSafe™ TURBO (40% in TURBO mode)	(9.2 lpm)

Made in USA


Product Ordering Information

Catalog Number:

PP012SEO 12 Per Carton

PlumePort™ SEO Laparoscopic Plume Filtration Device

Protecting healthcare workers and patients from the hazards of surgical smoke plume is particularly challenging when performing laparoscopic procedures. As plume begins to build within the peritoneal cavity, visibility concerns add to the safety hazards associated with surgical smoke plume. PlumePort SEO is a laparoscopic plume filtration device that utilizes available surgical plume evacuators and easily connects to standard trocars to safely remove the surgical plume during laparoscopy.



Buffalo Filter®

5900 Genesee Street
Lancaster, NY 14086
U.S.A.

Phone: 800.343.2324 (US)
Phone: +1.716.835.7000
Fax: 716.835.3414

Email: smoke-info@buffalofilter.com
Web: www.buffalofilter.com

PN 903280 REV B

Buffalo Filter® is a registered trademark of Buffalo Filter, LLC in the United States and other countries. All other trademarks used herein are the property of their respective owners. ©2013. Buffalo Filter. All rights reserved.