

# PlumeSafe® Turbo™

Surgical Smoke Plume Evacuator



## Interactive LCD Display

Our interactive display puts information at your fingertips.

## RFID Technology

Automatically tracks and records filter life.

## Turbo Mode

Significantly increases the airflow to quickly remove plume from the surgical site.

## Whisper Technology

Advanced sound control design ensures ultra-quiet operation with unprecedented sound quality.

## Simple Filter Changes



## Technical Information

### Noise Levels\*

@ maximum power < 58 dBA

### Dimensions (H x W x D)

13.5 in x 8 in x 17 in

34 cm x 20 cm x 43 cm

### Flow Rates\*

68 cfm (1923 lpm) using 1-1/4 in (32 mm) tubing

44 cfm (1257 lpm) using 7/8 in (22 mm) tubing

### Weight

With filter installed: 16.9 lbs (7.0 kg)

Without filter installed: 14.0 lbs (6.0 kg)

### Warranty

Three (3) Year Full Warranty

## Product Ordering Information

### Catalog Number:

PSWTURBO AS 100/120 Volt 50/60 Hz

PSWTURBO AS/220 220/240 Volt 50/60 Hz

(includes one footswitch; filter sold separately)

# PlumeSafe® Turbo™

## Surgical Smoke Plume Evacuator

PlumeSafe® TURBO™ is our most powerful smoke plume evacuator and features enhanced Whisper Technology, TURBO suction power, and an LCD display which provides a clear visual readout regarding exact filter time remaining, suction rate, mode, and filter type.

PlumeSafe® TURBO™ effectively captures and filters surgical smoke plume and removes odor and other potentially hazardous by-products. Each filter features RFID Technology and contains four stages of filtration in a single housing with a built-in pre-filter, special blend of activated carbon, ULPA filter and post filter. Operate with a foot pedal, control panel, or an automatic activation device.



Interactive LCD screen



### Buffalo Filter, LLC

5900 Genesee Street

Lancaster, NY 14086

U.S.A.

Phone: 800.343.2324 (US)

Phone: +1.716.835.7000

Fax: 716.835.3414

Email: [smoke-info@buffalofilter.com](mailto:smoke-info@buffalofilter.com)

Web: [www.buffalofilter.com](http://www.buffalofilter.com)

\*Nominal values

Made in USA

PN 903211 REV B



Pneumatic footswitch (included)



EZLink (sold separately; part: EZLINK01).

Buffalo Filter® is a registered trademark of Buffalo Filter, LLC, in the United States and other countries. All other trademarks used herein are the property of their respective owners. ©2014, Buffalo Filter. All rights reserved.