

Ø 2.3
mm

134 °C
273 °F

radial Spray Catheter

- for marking suspicious mucosa
- 360° radial spraying
- LL connector with mandrin

length	Order No.
120 cm	47 226-A
160 cm	47 226-B
180 cm	47 226-C
230 cm	47 226-G



Ø 1.8
mm

134 °C
273 °F

Spray Catheter for wide spray

- for marking suspicious mucosa
- for straight, cone-shaped spraying
- LL connector with mandrin

length	Order No.
120 cm	47 185-A
160 cm	47 185-B
180 cm	47 185-C
230 cm	47 185-G

