

ZepHr®

Impedance/pH Reflux Monitoring System

HOW ARE YOU ASSESSING YOUR PPI RESISTANCE PATIENTS TODAY?

Patients on PPI's who continue to experience symptoms such as cough, heartburn, regurgitation and chest pain often are difficult to diagnose using traditional acid (pH) monitoring approaches. In fact, a recent study¹ states that physicians using only acid (pH) monitoring for diagnostics, lack the capability of accurately diagnosing GERD in 35% of their patients.


The ZepHr® Impedance/pH Reflux Monitoring System employs impedance to detect ALL reflux activity and uses pH to categorize each episode as acid or nonacid for **Total Reflux Monitoring**. Comprehensive analysis quantifies all reflux patterns and symptom associations in patients studied on or off acid suppression medication.

Having introduced impedance/pH monitoring to the G.I. market, Sandhill continues to evolve this unprecedented technology... delivering all the information you need for a precise, comprehensive assessment of acid and nonacid reflux as well as the correlation between reflux and symptoms.

ALL THE DATA YOU NEED FOR PRECISE, ACCURATE ASSESSMENT OF YOUR PATIENTS.



ARE YOU CURRENTLY GETTING ALL THE DATA?



GET TOTAL REFLUX MONITORING ANALYSIS



SANDHILL UNIVERSITY - A HIGHER LEVEL OF TRAINING AND EDUCATION

Providing great products is just the first step in a long term relationship with our customer... we believe that education and training is the key component to the best utilization of our products which enables you to provide outstanding patient care. That's why at Sandhill University, we have developed a comprehensive set of options to meet your needs, including:

- Sandhill University – Denver
- Sandhill University – Online
- Webinars
- Webcasts
- Cybercoaching

Sandhill Scientific, Inc
9150 Commerce Center Circle, #500
Highlands Ranch, CO 80129 USA
800.468.4556 303.470.7020
Fax: 303.470.2975

Sandhill Scientific - UK Ltd.
Unit 13
Oasis Business Park
Eynsham
Oxfordshire
OX29 4TP
Tel: (+44) 1865-882532
Fax: (+44) 1865-883970

Sandhill Scientific - Europe
Decopra s.r.o
Svabky 2, Budova B
180 00 Praha 8
Czech Republic
Tel: (+420) 222 352 340
Fax: (+420) 266 310 110

* patent pending

1. Agrawal and Castell; Reflux ReductionTherapy: When PPI's Are Not Enough, Practical Gastroenterology, 68-72, April 2006.
2. Inder Mainie, R. Tutuiian & D. O. Castell; Medical University of South Carolina; Symptoms on PPI Therapy Associated with Nonacid, Acid or No Reflux; ACG Presentation; October 2004.
3. Patients with typical symptoms* (heartburn, regurgitation, chest pain) while on PPI Therapy.
US Patent #7,493,158

INDICATIONS FOR COMBINED IMPEDANCE/pH TESTING

- Persistent symptoms while on acid suppressive therapy
- Primarily postprandial symptoms
- Reflux symptoms and frequent meal ingestion (i.e., infant)



SMALL SIZE... BIG PERFORMANCE

Small by powerful! Your patients will appreciate the large, easy to understand controls including our well known symptom buttons that make reporting as easy as 1-2-3.

TREATMENT CONCLUSIONS


- Patients with nonacid identified by impedance/ pH whose symptoms have not responded to PPI therapy may benefit from the use of other medications.¹
- Clinical trials have established that nonacid reflux can be associated with GERD symptoms. In addition, ZepHr provides a true negative study by identifying patients with no reflux association.²
- Positive symptom index for nonacid or acid reflux using impedance/pH predicts successful response to Laparoscopic Nissen Fundoplication.²

ZepHr®

Impedance/pH Reflux Monitoring System

FEATURING ENHANCED BIOVIEW® G.I. DIAGNOSTIC SOFTWARE

Powered by Sandhill's revolutionary AutoSCAN®, BioVIEW analysis software simplifies interpretation of an impedance/pH study. Plus, BioVIEW easily guides you through a Z Swallow Challenge™ - a bolus transit screening test during an impedance/pH test with a single intubation that provides quantitative analysis of bolus transit effectiveness.



SyncVIEW™
uses
patented*
playback
to visualize
bolus transit

RapidVIEW™
quick
visualization
of an acid
or nonacid
reflux event

