

SURESTORE®

Endoscope Storage & Transportation




MEDICAL
INNOVATIONS
GROUP®

Medical Innovations Group are proud to introduce **SureStore®**, the Future of Endoscope Storage and Transportation.

British designed and manufactured by our experienced engineers and technical staff, **SureStore®** offers cutting edge long term "active" aseptic storage and transportation of Endoscopes (including GI and Theatre based scopes). Conditioned directly from the Endoscope Washer Disinfector (EWD) and completely negating the need for time consuming air drying, **SureStore®** validates each scope for storage of **up to 100 days!**

The unique design technology of **SureStore®** provides a 'Patient Ready' Packaged endoscope in Just 5 minutes.



SureStore® has been independently validated by the renowned Biotech Germande facility in accordance with the requirements of NFS 98-030: 2012 (prEN16442).

The System has been extensively tested with Endoscopes from Leading Manufacturers including Olympus, Pentax, Fujinon and Storz and results have shown no damaging or detrimental effect on the Endoscope.

KEY FEATURES INCLUDE:

- ✓ Storage Validation 100 Days
- ✓ IMS (Independent Monitoring System)
- ✓ Reduced Scope Wear & Tear
- ✓ Flexible Transportation
- ✓ Fully Traceable Cycle Data
- ✓ Quick & Simple
- ✓ Secure Robust Storage
- ✓ Space Saving



Space Saving



Quick & Simple



Fully Traceable

SureStore® enables scopes to be easily prepared in advance for emergency situations and also for any other out of hours, weekend/off-site endoscopy work, allowing for an amount of scopes to be made ready very quickly. The scopes can then be stored until needed, transported to various departments in the Hospital or to other clinics using the single **TravelScope®** case or the **TranScope** lockable cart.















The user friendly **SureStore®** conditioning unit requires no costly installation or endoscope drying facilities and is universal to any EWD. The purpose designed contoured tray and triple sealed vacuum pouch provide complete transport and storage protection for all types of endoscope.

SureStore® has the potential to significantly reduce staff time and endoscope wear and tear associated with reprocessing and transit.

TECHNICAL SPECIFICATIONS

- Process time - approximately 5 minutes from start to finish
- Dimensions
 - Height: 1220mm (960mm exc screen)
 - Width: 545mm
 - Depth: 690mm
- Installation space - unit does not require any additional ventilation space and can therefore be positioned against a wall or between other equipment.
- Validation - Biotech Germande and NFS 98-030
- Installation & Operating Requirements
 - Power Requirements 13amp 220 - 240V
 - Completely mobile on 4 independently lockable castors

ORDERING INFORMATION

	Product Ref	NHS SC*	Product Description	Units Per Box
	100600	FAL13614	SureStore® Set Includes: 1 x SureStore® Conditioning Unit 12 x SureStore® Supportive Trays 1 x 5 Litre N-Sure® Solution 50 x SureStore® Outer Pouches 45 x SureStore® Sterile Collection Pads 50 x Cleanascope® Liners	1
	100602	FAL13615	N-Sure® Solution for use with SureStore® endoscope conditioning system 5L	1
	100603	FAL13616	N-Sure® Solution for use with SureStore® endoscope conditioning system 5L	4
	100604	FAL13617	SureStore Outer Pouch	100
	100605	FAL13629	SureStore® Sterile Collection Pads	225
	100606	FAL13620	SureStore® Supportive Tray The SureStore® Supportive Tray has been specifically contoured for the protection of conditioned flexible endoscopes *Fully compatible with the SureStore® Conditioning System	1
	100607	FAL13621	SureStore® Supportive Tray	5
	100608	FAL13622	SureStore® Supportive Tray	10
	700344	FAL13618	SureStore® Tracking Labels for use with SureStore® conditioning unit	400
	200018	FAL13623	SureStore® Hepa Air Filter	1
	103580	FAL12799	CleanaScope® sterile liner set (red and green covers with base liner)	200
	100521	FAL13619	TravelScope® Transit Case - The TravelScope® transit case offers complete protection and flexibility in transporting fully conditioned flexible endoscopes within the healthcare facility and beyond	1
	103543	-	TranScope Lockable Cart 6 Tray	1
	103544	-	TranScope Lockable Cart 6 Tray with Lids	1
	105047	-	Virusolve+ Eco Wipes 225 Sheets	4
	105045	-	Virusolve+ Trigger Spray 750ml	12

*Available via NHS Supply Chain (NHS SC)

+44 (0)1702 291878
 Surestore@medinn.co.uk
 sure-store.co.uk
 Campfield Road, Shoeburyness
 Essex, SS3 9BX, UK

