

PlumePen® Pro

Surgical Smoke Evacuation Pencil



BUFFALO FILTER®

Filtration Group®

Sleek, Slim Profile

Smaller than competitive surgical smoke evacuation pencils; similar in size to standard ESU pencils

Removable Smoke Capture Port

Allows for easy electrode and capture port replacement

Compatible with Standard Electrodes

Designed for surgeries where longer electrodes are desired

Full Blade Tip Reveal

Provides an unobstructed view of surgical site

360° Swivel

Allows for fluid hand movement

Textured Profile

Improves surgeon's grip



Actual Size

Additional capture ports available to accommodate 4 in (102 mm) and 6 in (152 mm) electrodes.

Product Ordering Information

Pencils:

PlumePen Pro (Non-Stick ESU Blade)
PLPRO4020 10 ft (3 m) tubing 20/Case
PLPRO4520 15 ft (4.6 m) tubing 20/Case

Capture Ports (X-ray detectable):

PROTIP420 20/Case
(Accommodates 4 in [102 mm] electrodes)
PROTIP620 20/Case
(Accommodates 6 in [152 mm] electrodes)

PlumePen® Pro

Surgical Smoke Evacuation Pencil

The new, improved PlumePen Pro features an easy to remove x-ray detectable capture port that allows for effortless electrode replacement to accommodate surgeries where longer electrodes are desired. Pro is smaller than competitive surgical smoke evacuation pencils and similar in size to standard electrosurgical pencils.

Pro is a cost effective and suitable alternative for surgeries where a pencil adapter is less optimal.

**Capture
More Smoke**
Up to 32% higher flow rate
than competitive pencils.



Additional capture port lengths available depending on the type of surgery and length of electrode used.



(each sold separately)



Buffalo Filter LLC

5900 Genesee Street
Lancaster, NY 14086
U.S.A.

Phone: 800.343.2324 (US)
Phone: +1.716.835.7000
Fax: 716.835.3414

Email: smoke-info@buffalofilter.com
Web: www.buffalofilter.com