



STERILE EO

### Biliary Single-Size Dilation Balloon

- total length 200 cm
- radiopaque markers



Order No.	Ø balloon	balloon length	fitting to working channel	Ø catheter	accepts guide wire
341 0420	4.0 mm	20.0 mm	2.8 mm	7 FR	.035 inch
341 0620	6.0 mm	20.0 mm	2.8 mm	7 FR	.035 inch
341 0830	8.0 mm	30.0 mm	2.8 mm	7 FR	.035 inch
341 1030	10.0 mm	30.0 mm	2.8 mm	7 FR	.035 inch
341 1230	12.0 mm	30.0 mm	2.8 mm	7 FR	.035 inch

### Esophageal Single-Size Dilation Balloon

- total length 200 cm
- radiopaque markers



Order No.	Ø balloon	balloon length	fitting to working channel	Ø catheter	accepts guide wire
341 0880	8.0 mm	80.0 mm	2.8 mm	7 FR	.035 inch
341 1280	12.0 mm	80.0 mm	2.8 mm	7 FR	.035 inch
341 1680	16.0 mm	80.0 mm	3.2 mm	7 FR	.035 inch
341 1880	18.0 mm	80.0 mm	3.2 mm	7 FR	.035 inch
340 0880	8.0 mm	80.0 mm	2.8 mm	7 FR	
340 1280	12.0 mm	80.0 mm	2.8 mm	7 FR	
340 1680	16.0 mm	80.0 mm	3.2 mm	7 FR	
340 1880	18.0 mm	80.0 mm	3.2 mm	7 FR	

### Achalasia Single-Size Balloon Dilator

- total length 100 cm



Order No.	Ø balloon	balloon length	fitting to working channel	Ø catheter	accepts guide wire
341 3080	30.0 mm	80.0 mm		4.7 mm	.035 inch
341 3580	35.0 mm	80.0 mm		4.7 mm	.035 inch
341 4080	40.0 mm	80.0 mm		4.7 mm	.035 inch

## recommended equipment

Inflation / Deflation Device  
incl. PSI-Manometer 0-440 PSI



STERILE EO

Order No.

340 800 INFL

