


single-use

self-expanding Biliary Stents



### Characteristics

- TTS (Through-The-Scope)
- Application Set
  - non-coated stents Fr. 8
  - coated stents Fr. 9
- Length: 180 cm
- Recommended guide wire: .035 inch (0.9 mm)

### Placement method

- Proximal

### Special features

- braided Teflon tubing with transparent distal end
- precise trans-papillary Stent placement under constant endoscopic view
- no shortening at proximal end of the stent (close to the sphincter)
- minimum / less fluoroscopy needed
- easy handling and optimal deployment time
- capacity for repositioning of up to 60%
- indicator for "point of no return"



STERILE EO

### self-expanding Biliary Stents

- incl. pre-loaded application set
- trans-papillary (proximal) placement
- radiopaque markers

#### uncoated



Order No.	length	stent diameter
BIL-0-10- 40-RP	40 mm	10 mm
BIL-0-10- 60-RP	60 mm	10 mm
BIL-0-10- 80-RP	80 mm	10 mm
BIL-0-10-100-RP	100 mm	10 mm

#### silicone coated



Order No.	length	stent diameter
BIL-1-10-40-RP	40 mm	10 mm
BIL-1-10-60-RP	60 mm	10 mm
BIL-1-10-80-RP	80 mm	10 mm