

### Characteristics

- OTW (Over-The-Wire)
- Application Set 4 or 6 mm, depending on stent diameter
- Length: 60 cm
- Recommended guide wire: .035 inch (0.9 mm)

### Placement method

- Distal

### Special features


- transparent outer tubing with graduation
- easy handling and optimal deployment time
- stent design with smooth atraumatic stent ends
- effective anchorage and migration prevention



STERILE EO

### self-expanding bronchial / tracheal Stents

- incl. pre-loaded application set
- distal placement
- radiopaque markers




### uncoated

Order No.	length	stent diameter
BRT-0-12-40	40 mm	12 mm
BRT-0-18-60	60 mm	18 mm

### self-expanding bronchial / tracheal Stents

- incl. pre-loaded application set
- distal placement
- radiopaque markers



### silicone coated

Order No.	length	stent diameter
BRT-2-10-30	30 mm	10 mm
BRT-2-20-70	70 mm	20 mm

Due to the variety of possibly needed lengths and diameters, bronchial/tracheal stents are mainly produced on demand. The above codes are meant as ordering examples, for other available dimensions please contact ENDO-FLEX or your local distribution partner.