




Multi-Stage-Balloon, Latex-Free

- 3 Diameter
- for guide wire .035 inch
- improved radiopaque markers
- three different syringes included



triple lumen



catheter length	catheter	Ø balloon inflated (mm)			Order No.
		1.5 ml	2.5 ml	5.0 ml	
200 cm	7 Fr.	9 mm	13 mm	16 mm	N135030PRO

double lumen



catheter length	catheter	Ø balloon inflated (mm)			Order No.
		1.5 ml	2.5 ml	5.0 ml	
200 cm	7 Fr.	9 mm	13 mm	16 mm	N135010PRO



single-use

Extraction Balloons