


polypectomy snare



An extensive line of polypectomy snares...

From cold snaring to advanced polypectomy, US Endoscopy offers a product solution to fit a variety of clinical needs.

product information

product number	description	needle gauge/ projection	sheath diameter (mm)	length (cm)	approximate snare size (cm)	units/box
00711018	Short Throw™ snare - mini	-	2.4	230	1.25 x 3.0	10
00711019	Short Throw™ snare - standard oval	-	2.4	230	2.5 x 5.4	10
00711120	Rotator™ snare - standard	-	2.4	230	2.5 x 5.4	10
00711121	Rotator™ snare - mini	-	2.4	230	1.25 x 3.0	10

All active cord connections are Olympus style

For more information on the polypectomy snares, contact US Endoscopy today.

Distributed by:



listening...and delivering solutions

5976 Hefley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@userendoscopy.com

www.userendoscopy.com

760908 Rev. B

US Endoscopy is a subsidiary of STERIS Corporation.
Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™
are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

Follow us:    