




**Fr. 7** → **Fr. 5**  Triple lumen tube distally tapered from 2.2 mm to 1.8 mm in treatment area.


**triple lumen** 

**Tapered Tip**

- for guide wire .035 inch




	braided cutting wire			monofilament cutting wire		
packing unit	20 mm	25 mm	30 mm	20 mm	25 mm	30 mm
5 pieces	OE 104 2220TRL	OE 104 2225TRL	OE 104 2230TRL	OE 104 2220TRL-M	OE 104 2225TRL-M	OE 104 2230TRL-M


**triple lumen** 

**Standard Tip**

- for guide wire .035 inch




	braided cutting wire			monofilament cutting wire		
packing unit	20 mm	25 mm	30 mm	20 mm	25 mm	30 mm
5 pieces	OE 100 2220TRL	OE 100 2225TRL	OE 100 2230TRL	OE 100 2220TRL-M	OE 100 2225TRL-M	OE 100 2230TRL-M


**triple lumen** 

**Smart-Cut**

- partially insulated cutting wire
- for guide wire .035 inch




packing unit	20 mm	25 mm	30 mm
5 pieces	OE 104 2220TRLC	OE 104 2225TRLC	OE 104 2230TRLC

**triple lumen** 

**Precut-Type**

- for guide wire .035 inch



packing unit	20 mm	25 mm	30 mm
5 pieces	OE 102 2220TRL	OE 102 2225TRL	OE 102 2230TRL



**Fr. 7**



**Fr. 5**



STERILE EO

Double lumen tube distally tapered from 2.2 mm to 1.8 mm in treatment area.

double lumen



Tapered Tip

- for guide wire .035 inch



packing unit	20 mm	25 mm	30 mm
5 pieces	OE 104 2220DL	OE 104 2225DL	OE 104 2230DL

double lumen



Smart-Cut

- partially insulated cutting wire
- for guide wire .035 inch



packing unit	20 mm	25 mm	30 mm
5 pieces	OE 104 2220DLC	OE 104 2225DLC	OE 104 2230DLC

double lumen



Standard Tip

- for guide wire .035 inch



packing unit	20 mm	25 mm	30 mm
5 pieces	OE 100 2220DL	OE 100 2225DL	OE 100 2230DL

double lumen



Precut-Type

- for guide wire .035 inch



packing unit	20 mm	25 mm	30 mm
5 pieces	OE 102 2220DL	OE 102 2225DL	OE 102 2230DL



Fr. 5



STERILE EO

single lumen



Tapered Tip

- no guide wire access



packing unit	20 mm	30 mm
5 pieces	OE 105 1821	OE 105 1831

single lumen



Precut-Type

- no guide wire access



packing unit	20 mm	30 mm
5 pieces	OE 103 1821	OE 103 1831

single lumen



Precut-Type with Tapered Tip

- no guide wire access



packing unit	20 mm	30 mm
5 pieces	OE 105 1821PC	OE 105 1831PC



Fr. 7



Fr. 5



Double and Triple lumen tubes distally tapered from 2.2 mm to 1.8 mm in treatment area.

triple lumen



Needle-Knife, Pointed-Type

- for guide wire .035 inch

stop at 6 mm

packing unit

5 pieces	OE 110 18N-TRL
----------	----------------



double lumen



Needle-Knife, Pointed-Type

- for guide wire .035 inch

stop at 6 mm

packing unit

5 pieces	OE 110 18N-DL
----------	---------------



Fr. 5



single lumen



Needle-Knife, Pointed-Type

- no guide wire access

stop at 3 mm

packing unit

5 pieces	OE 110 18N3
----------	-------------



Needle-Knife, Blunt-Type

- no guide wire access

stop at 3 mm

packing unit

5 pieces	OE 110 18SP3
----------	--------------



**Fr. 7**

STERILE EO

single lumen


B-II-Type

- for guide wire .035 inch

packing unit

5 pieces

OE 112 22FD

**Fr. 5**

STERILE EO

single lumen


B-II-Type

- no guide wire access

packing unit

5 pieces

OE 112 18



single-use

B-II-Type Sphincterotomes