

# Flexible Endoscopy CartOPF-1



Twin Scope Holder  
(optional)



Range of Monitor Posts or  
Articulating Arms available  
(optional)



Sterile Water  
Bottle Holder  
(optional)




Keyboard Drawer  
W: 460mm x D: 230mm  
x H: 60mm  
(standard)



Pushing Handle  
(standard)



Multi Size Gas  
Bottle Holder  
(optional)



# Flexible Endoscopy Cart OPF-1

The Flexible Endoscopy Cart OPF-1 has been designed to suit all budgets and is compatible with most Olympus, Pentax and Fuji equipment.

Supplied as standard with 2 adjustable equipment shelves, 1 keyboard drawer, a large pushing handle, ABS top tray, 6 way socket outlet and 4 x 5" castors with full buffering system (2 lockable). The cart can be fully customised with any of the optional accessory items listed below.

**Cart weight:** 42kg

## Accessory Items

